

Anti-Nyx Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Gzu5

Item Number: CM0025472

Introduction

Rabbit polyclonal antibody targeting human Nyctalopin (NYX), reactive with human and mouse. Suitable for Western blot, immunohistochemistry (paraffin), and ELISA. Ideal for retinal biology and vision disorder studies.

[Learn More](#)

Parameter	Value
Product Name	NYX Rabbit pAb
Remarks/Alias	CLRP; NBM1; CSNB1; CSNB4; CSNB1A; NYX
Species	Human
GeneID (Human)	60506
GeneID	60506
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 282-481 of human NYX (NP_072089.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9GZU5
Protein Weight	52kDa
Shipping	Ice bag