

Anti-Oaf Polyclonal Antibody For Wb, Elisa - Q86Ud1

Item Number: CM0014758

Introduction

Rabbit polyclonal antibody targeting human OAF, validated for Western blot and ELISA. Cross-reacts with mouse and rat OAF. Suitable for viral research and proteomics applications.

[Learn More](#)

Parameter	Value
Product Name	OAF Rabbit pAb
Remarks/Alias	NS5ATP13TP2; OAF
Species	Human
GeneID (Human)	220323
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-150 of human OAF (NP_848602.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q86UD1
Protein Weight	31kDa
Shipping	Ice bag