

Anti-Oas1 Polyclonal Antibody For Wb, Elisa - P00973

Item Number: CM0020953



Introduction

Rabbit polyclonal antibody targeting human OAS1 (2'-5'-oligoadenylate synthase 1), validated for Western blot and ELISA. Reacts with human and mouse. Essential for innate antiviral immunity research and IVD development.

[Learn More](#)

Parameter	Specification
Product Name	OAS1 Rabbit pAb
Remarks/Alias	OIAS; IFI-4; OIASI; IMD100; E18/E16; OAS1
Species	Human
GeneID	4938
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-364 of human OAS1 (NP_058132.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00973
Protein Weight	46kDa
Shipping	Ice bag