

Anti-Oas1 Rabbit Polyclonal Antibody For Wb And Elisa - P00973

Item Number: CM0015205



Introduction

OAS1 rabbit polyclonal antibody validated for Western blot (WB) and ELISA applications. Detects human OAS1 with cross-reactivity to mouse and rat. Ideal for antiviral innate immunity and 2-5A pathway research.

[Learn More](#)

Parameter	Value
Product Name	OAS1 Rabbit pAb
Remarks/Alias	OIAS; IFI-4; OIASI; IMD100; E18/E16; OAS1
Species	Human
Gene ID (Human)	4938
Gene ID	4938
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human OAS1 (NP_058132.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P00973
Protein Weight	46kDa
Shipping	Ice bag