

Anti-Oas3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6K5

Item Number: CM0026785



Introduction

OAS3 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Cross-reactive with human, mouse, and rat. Suitable for innate immunity and antiviral research. 121 kDa protein weight.

[Learn More](#)

Parameter	Value
Product Name	OAS3 Rabbit pAb
Remarks/Alias	p100; p100OAS; OAS3
Species	Human
GeneID (Human)	4940
GeneID	4940
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-620 of human OAS3 (NP_006178.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6K5
Protein Weight	121kDa
Shipping	Ice bag