


Anti-Oasl Rabbit Polyclonal Antibody For Wb, Elisa - Q15646

Item Number: CM0022738



Introduction

Rabbit polyclonal antibody against human OASL (2'-5'-oligoadenylate synthetase-like protein), validated for Western blot and ELISA. Cross-reacts with human and mouse OASL. Ideal for antiviral immunity research.

[Learn More](#)

Parameter	Specification
Product Name	OASL Rabbit pAb
Remarks/Alias	OASL; OASLd; TRIP-14; TRIP14; p59 OASL; p59-OASL; p59OASL; 2'-5'-oligoadenylate synthetase like
Species	Human
Gene ID	8638
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human OASL (NP_003724.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15646
Protein Weight	29kDa/43kDa/59kDa
Shipping	Ice bag