

Anti-Oat Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P04181

Item Number: CM0009154



Introduction

Rabbit monoclonal antibody against human OAT (Ornithine aminotransferase), validated for WB, IF, IHC, ELISA. Reacts with human, mouse, rat. Suitable for studying amino acid metabolism and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] OAT Rabbit mAb
Remarks/Alias	OKT; GACR; HOGA; OATASE
Species	Human
Gene ID (Human)	4942
Gene ID	4942
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 160-439 of human OAT (NP_000265.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04181
Protein Weight	49kDa
Shipping	Ice bag