

Anti-Oat Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P04181

Item Number: CM0026423



Introduction

Rabbit polyclonal antibody targeting human OAT (SWISS P04181). Validated for WB, IHC-P, IP, and ELISA. Detects OAT in human, mouse, and rat. Suitable for studying amino acid metabolism and mitochondrial diseases.

[Learn More](#)

Parameter	Specification
Product Name	OAT Rabbit pAb
Remarks/Alias	OKT; GACR; HOGA; OATASE; OAT
Species	Human
Gene ID (Human)	4942
Gene ID	4942
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04181
Protein Weight	49kDa
Shipping	Ice bag