

Anti-Oat3 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q8Tcc7

Item Number: CM0024439

Introduction

Rabbit polyclonal antibody against human OAT3, validated for IF/ICC and ELISA. Detects the basolateral organic anion transporter SLC22A8, crucial for renal drug excretion and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	OAT3 Rabbit pAb
Remarks/Alias	OAT3
Species	Human
Gene ID (Human)	9376
Gene ID	9376
Immunogen	Recombinant protein corresponding to a sequence within amino acids 32-123 of human OAT3 (NP_004245.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q8TCC7
Protein Weight	60kDa
Shipping	Ice bag