

Anti-Oat3/Slc22A8 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tcc7

Item Number: CM0025497

Introduction

This rabbit polyclonal antibody specifically detects human OAT3/SLC22A8 protein. Validated for Western blot (WB) and ELISA, it targets amino acids 443-542. Ideal for renal transporter and drug disposition research.

[Learn More](#)

Parameter	Value
Product Name	OAT3/SLC22A8 Rabbit pAb
Remarks/Alias	OAT3; OAT3/SLC22A8
Species	Human
GeneID (Human)	9376
GeneID	9376
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 443-542 of human OAT3/SLC22A8 (NP_004245.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TCC7
Protein Weight	60kDa
Shipping	Ice bag