

Anti-Oaz1 Rabbit Polyclonal Antibody For Wb, Elisa - P54368

Item Number: CM0029780



Introduction

OAZ1 Rabbit Polyclonal Antibody validated for WB and ELISA, reactive with Human, Mouse, Rat. Regulates polyamine metabolism and ODC degradation. Recombinant immunogen. ~25 kDa.

[Learn More](#)

Parameter	Value
Product Name	OAZ1 Rabbit pAb
Remarks/Alias	AZ1; AZI; OAZ; OAZ1
Species	Human
Gene ID	4946
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human OAZ1 (NP_004143.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54368
Protein Weight	~25 kDa
Shipping	Ice bag