


Anti-Obecn Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q5Vst9

Item Number: CM0021061



Introduction

Rabbit polyclonal antibody targeting human OBSCN, applicable for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for striated muscle and sarcomere research.

[Learn More](#)

Parameter	Value
Product Name	OBSCN Rabbit pAb
Remarks/Alias	UNC89; RHABDO1; ARHGEF30; OBSCN
Species	Human
GeneID (human)	84033
GeneID	84033
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1621-1710 of human OBSCN (NP_443075.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5VST9
Protein Weight	868kDa
Shipping	Ice bag