

Anti-Oca2 Rabbit Polyclonal Antibody For Wb, Elisa - Q04671

Item Number: CM0027603

Introduction

Rabbit polyclonal antibody against human OCA2/P protein, validated for WB and ELISA. Cross-reacts with mouse. Suitable for melanogenesis, pigmentation, and eye color research. Shipped on ice.

[Learn More](#)

Parameter	Value
Product Name	OCA2 Rabbit pAb
Remarks/Alias	P, BEY, PED, BEY1, BEY2, BOCA, EYCL, HCL3, EYCL2, EYCL3, SHEP1, D15S12, OCA2
Species	Human
GeneID (Human)	4948
GeneID	4948
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human OCA2 (NP_000266.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q04671
Protein Weight	93kDa
Shipping	Ice bag