

Anti-Occludin Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q16625

Item Number: CM0030566



Introduction

Occludin rabbit polyclonal antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA with cross-reactivity to human, mouse, and rat tight junction protein for paracellular permeability and HCV entry research.

[Learn More](#)

Parameter	Value
Product Name	Occludin Rabbit pAb
Remarks/Alias	BLCPMG; PTORCH1; PPP1R115; Occludin
Species	Human
GeneID (Human)	100506658
GeneID	100506658
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-522 of human Occludin (NP_002529.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16625
Protein Weight	59kDa
Shipping	Ice bag