

Anti-Ociad1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nx40

Item Number: CM0017808

Introduction

Rabbit polyclonal antibody against human OCIAD1, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for stem cell and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	OCIAD1 Rabbit pAb
Remarks/Alias	OCIA; ASRIJ; TPA018; OCIAD1
Species	Human
GeneID (Human)	54940
GeneID	54940
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-245 of human OCIAD1 (NP_060300.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NX40
Protein Weight	28kDa
Shipping	Ice bag