

Anti-Ociad2 Polyclonal Antibody For Wb - Q56VI3

Item Number: CM0015620

Introduction

High-quality OCIAD2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human samples. Ideal for mitochondrial complex III assembly and STAT3 signaling research.

[Learn More](#)

Specification	Value
Product Name	OCIAD2 Rabbit pAb
Remarks/Alias	OCIAD2
Species	Human
GeneID	132299
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-99 of human OCIAD2 (NP_689611.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q56VL3
Protein Weight	17kDa
Shipping	Ice bag