

Anti-Ocrl Rabbit Polyclonal Antibody For Wb And Elisa - Q01968

Item Number: CM0022459



Introduction

Rabbit polyclonal antibody against human OCRL (Q01968), validated for WB and ELISA. Detects human, mouse, and rat OCRL; ideal for studying inositol phosphate metabolism, endosomal trafficking, and ciliopathies.

[Learn More](#)

Parameter	Value
Product Name	OCRL Rabbit pAb
Remarks/Alias	LOCR; DENT2; NPHL2; OCRL1; Dent-2; INPP5F; OCRL-1; OCRL
Species	Human
Gene ID (Human)	4952
Gene ID	4952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human OCRL (NP_001578.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01968
Protein Weight	104kDa
Shipping	Ice bag