

Anti-Odc1 Rabbit Monoclonal Antibody For Wb, Elisa - P11926

Item Number: CM0008415



Introduction

Rabbit monoclonal antibody against human ODC1 (ornithine decarboxylase), validated for Western blotting (WB) and ELISA. Detect a band at ~51 kDa. Ideal for studying polyamine biosynthesis and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	ODC1 Rabbit mAb
Remarks/Alias	ODC; BABS; NEDBA; NEDBIA; ODC1
Species	Human
GeneID (Human)	4953
GeneID	4953
Immunogen	Synthetic Peptide. A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ODC1 (P11926).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P11926
Protein Weight	51kDa
Shipping	Ice bag