

## Anti-Odc1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P11926

Item Number: CM0019532



### Introduction

Rabbit polyclonal antibody against human ODC1 (Ornithine decarboxylase), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for polyamine biosynthesis and cell proliferation research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | ODC1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | ODC; BABS; NEDBA; NEDBIA; ODC1                                                                                  |
| Species          | Human                                                                                                           |
| Gene ID (Human)  | 4953                                                                                                            |
| Gene ID          | 4953                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ODC1 (NP_002530.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IF/ICC, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | P11926                                                                                                          |
| Protein Weight   | 51kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |