

Anti-Odf1 Polyclonal Antibody For Wb And Elisa - Q14990

Item Number: CM0022713

Introduction

Rabbit polyclonal antibody against human ODF1 for WB and ELISA, cross-reactive with mouse and rat, suitable for sperm tail outer dense fiber and spermatogenesis research, protein weight 28kDa.

[Learn More](#)

Parameter	Value
Product Name	ODF1 Rabbit pAb
Remarks/Alias	ODF; RT7; ODF2; ODFP; SODF; CT133; ODF27; ODFPG; HSPB10; ODFPGA; ODFPGB; ODF1
Species	Human
GeneID (Human)	4956
GeneID	4956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-180 of human ODF1 (NP_077721.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14990
Protein Weight	28kDa
Shipping	Ice bag