

Anti-Odf2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Bjf6

Item Number: CM0016524

Introduction

Rabbit polyclonal antibody against human ODF2, validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat ODF2 (~95 kDa). Suitable for centrosome/centriole and sperm tail research.

[Learn More](#)

Parameter	Value
Product Name	ODF2 Rabbit pAb
Remarks/Alias	CT134; ODF84; ODF2/1; ODF2/2; ODF2
Species	Human
GeneID (Human)	4957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 630-829 of human ODF2 (NP_702911.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5BJF6
Protein Weight	95kDa
Shipping	Ice bag