

Anti-Ogdh Monoclonal Antibody For Wb, Ihc-P, Elisa - Q02218

Item Number: CM0003151



Introduction

Rabbit monoclonal antibody against human OGDH, suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studies on Krebs cycle and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	OGDH Rabbit mAb
Remarks/Alias	E1k; E1o; KGD1; OGDC; AKGDH; OGDH2; OGDHD; OGDH-E1; OGDH
Species	Human
GeneID (Human)	4967
GeneID	4967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 869-1023 of human OGDH (NP_002532.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02218
Protein Weight	116 kDa
Shipping	Ice bag