


Anti-Ogdh Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q02218

Item Number: CM0030231



Introduction

OGDH Rabbit pAb is a polyclonal antibody for detecting 2-oxoglutarate dehydrogenase E1 (OGDH) in human, mouse, rat samples via WB, IHC-P, IF/ICC, ELISA. Suitable for Krebs cycle and mitochondrial metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	OGDH Rabbit pAb
Remarks/Alias	E1k; E1o; KGD1; OGDC; AKGDH; OGDH2; OGDHD; OGDH-E1; OGDH
Species	Human
GeneID (Human)	4967
GeneID	4967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 148-427 of human OGDH (NP_001003941.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02218
Protein Weight	116kDa
Shipping	Ice bag