


## Anti-Ogdh Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q02218

Item Number: CM0016329



### Introduction

Anti-OGDH Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat OGDH, a key mitochondrial enzyme in the citric acid cycle. Suitable for metabolism and mitochondrial research.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | OGDH Rabbit pAb                                                                                                      |
| Remarks/Alias    | E1k; E1o; KGD1; OGDH; AKGDH; OGDH2; OGDHD; OGDH-E1; OGDH                                                             |
| Species          | Human                                                                                                                |
| Gene ID          | 4967                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 148-427 of human OGDH (NP_001003941.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | Q02218                                                                                                               |
| Protein Weight   | 116kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                              |