

Anti-Ogfr Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nzt2

Item Number: CM0025478

Introduction

Rabbit polyclonal antibody targeting human OGFR (opioid growth factor receptor). Suitable for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Designed for research on growth regulation and opioid signaling.

[Learn More](#)

Parameter	Value
Product Name	OGFR Rabbit pAb
Remarks/Alias	OGFR
Species	Human
GenelD (Human)	11054
GenelD	11054
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 358-482 of human OGFR (NP_031372.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZT2
Protein Weight	73kDa
Shipping	Ice bag