


Anti-Ogg1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O15527

Item Number: CM0008061



Introduction

OGG1 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA in human, rat. Targets human OGG1 (O15527), a DNA repair enzyme excising 8-oxoguanine lesions. Ideal for cancer and aging studies.

[Learn More](#)

Parameter	Value
Product Name	OGG1 Rabbit mAb
Remarks/Alias	HMMH; MUTM; OGH1; HOGG1; OGG1
Species	Human
GeneID (Human)	4968
GeneID	4968
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 246-345 of human OGG1 (O15527).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O15527
Protein Weight	39 kDa
Shipping	Ice bag