


Anti-Ogg1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O15527

Item Number: CM0017670



Introduction

Rabbit polyclonal antibody against human OGG1 (8-oxoguanine DNA glycosylase). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Use for DNA repair studies targeting the N-glycosylase/DNA lyase.

[Learn More](#)

Parameter	Value
Product Name	OGG1 Rabbit pAb
Remarks/Alias	HMMH; MUTM; OGH1; HOGG1; OGG1
Species	Human
Gene ID (Human)	4968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-345 of human OGG1 (NP_002533.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15527
Protein Weight	39kDa
Shipping	Ice bag