


Anti-Ogt Monoclonal Antibody For Wb, Elisa, Chip - O15294

Item Number: CM0005656



Introduction

Rabbit monoclonal antibody targeting O-GlcNAc transferase (OGT), validated for Western Blot, ELISA, and ChIP, with cross-reactivity to human, mouse, and rat. Ideal for studying O-GlcNAcylation, chromatin regulation, and insulin signaling.

[Learn More](#)

Parameter	Value
Product Name	OGT Rabbit mAb
Remarks/Alias	OGT1; HRNT1; MRX106; XLID106; HINCUT-1; O-GLCNAC; OGT
Species	Human
GeneID (human)	8473
GeneID	8473
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 947-1046 of human OGT (O15294).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	O15294
Protein Weight	117kDa
Shipping	Ice bag