


Anti-Ola1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Ntk5

Item Number: CM0016822



Introduction

Rabbit polyclonal antibody against human OLA1 (DOC45, GTPBP9). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for cancer research and cellular stress studies.

[Learn More](#)

Parameter	Value
Product Name	OLA1 Rabbit pAb
Remarks/Alias	DOC45, GBP45, GTPBP9, GTPBP9, PTD004, OLA1
Species	Human
GeneID	29789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-310 of human OLA1 (NP_037473.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NTK5
Protein Weight	45kDa
Shipping	Ice bag