


Anti-Olfm1 Rabbit Polyclonal Antibody For Wb And Elisa - Q99784

Item Number: CM0018060



Introduction

Rabbit polyclonal antibody against human OLFM1 (Noelin). Validated for WB and ELISA. Reacts with human and mouse. Useful for studying axonal growth regulation and neural development.

[Learn More](#)

Parameter	Value
Product Name	OLFM1 Rabbit pAb
Remarks/Alias	AMY; NOE1; OlfA; NOELIN1; OLFM1
Species	Human
Gene ID (Human)	10439
Gene ID	10439
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99784
Protein Weight	55kDa
Shipping	Ice bag