

Anti-Olfm3 Rabbit Polyclonal Antibody For Wb And Elisa - Olfm3 Rabbit Pab

Item Number: CM0019311

Introduction

Polyclonal anti-OLFM3 antibody raised in rabbit against a synthetic peptide from mouse OLFM3. Suitable for Western blot and ELISA detection of OLFM3 in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	OLFM3 Rabbit pAb
Remarks/Alias	B230206G02Rik; OLFM3
Species	Human
Gene ID (Human)	229759
Gene ID	229759
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 200-300 of mouse OLFM3 (NP_694797.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag