

Anti-Olfm4 Rabbit Polyclonal Antibody For Wb And Elisa - Q6Ux06

Item Number: CM0014623

Introduction

Rabbit polyclonal antibody against human OLFM4, validated for Western blotting and ELISA. Detects OLFM4 in human, mouse, and rat samples. Suitable for cancer research and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	OLFM4 Rabbit pAb
Remarks/Alias	GCL1; OLM4; OlfD; GW112; hGC-1; hOLF4; UNQ362; bA209J19.1; OLFM4
Species	Human
GeneID (Human)	10562
GeneID	10562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 211-510 of human OLFM4 (NP_006409.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UX06
Protein Weight	57kDa
Shipping	Ice bag