

Anti-Olig1 Monoclonal Antibody For Wb, Elisa - Q8Tak6

Item Number: CM0002173

Introduction

Rabbit monoclonal antibody targeting human Olig1 (Q8TAK6). Validated for Western blot and ELISA. Recognizes human, mouse, and rat Olig1. Suitable for research in gliogenesis, oligodendrocyte development, and glioma biology.

[Learn More](#)

Parameter	Value
Product Name	Olig1 Rabbit mAb
Remarks/Alias	BHLHB6; BHLHE21; Olig1
Species	Human
GeneID	116448
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 172-271 of human Olig1 (Q8TAK6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TAK6
Protein Weight	28kDa
Shipping	Ice bag