

Anti-Olig2 Monoclonal Antibody For Wb, If, Ihc, Elisa - Q13516

Item Number: CM0008064

Introduction

Olig2 Rabbit Monoclonal Antibody validated for WB, IF-P, IHC-P, ELISA. Detects endogenous levels of Olig2 in human, mouse, and rat. Suitable for oligodendrocyte, motor neuron, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Olig2 Rabbit mAb
Remarks/Alias	BHLHB1; OLIGO2; RACK17; PRKCBP2; bHLHe19
Species	Human
GeneID (human)	10215
GeneID	10215
Immunogen	Synthetic peptide of human Olig2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13516
Protein Weight	32 kDa
Shipping	Ice bag