

Anti-Olig2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13516

Item Number: CM0014733

Introduction

Rabbit polyclonal antibody targeting human Olig2, suitable for Western blot, immunohistochemistry, and ELISA. Detects human, mouse, and rat Olig2. Ideal for neurodevelopmental and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Olig2 Rabbit pAb
Remarks/Alias	BHLHB1, OLIGO2, RACK17, PRKCBP2, bHLHe19, Olig2
Species	Human
Gene ID (Human)	10215
Gene ID	10215
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-80 of human Olig2 (NP_005797.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13516
Protein Weight	32kDa
Shipping	Ice bag