

Anti-Olig2 Rabbit Monoclonal Antibody For Wb, Elisa - Q13516

Item Number: CM0010100

Introduction

Olig2 rabbit monoclonal antibody validated for Western blot and ELISA. Detects endogenous human, mouse, and rat Olig2 (~32 kDa). Essential for oligodendrocyte and motor neuron development research.

[Learn More](#)

Parameter	Value
Product Name	Olig2 Rabbit mAb
Remarks/Alias	BHLHB1; OLIGO2; RACK17; PRKCBP2; bHLHe19; Olig2
Species	Human
GeneID (Human)	10215
GeneID	10215
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Olig2 (Q13516).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13516
Protein Weight	32kDa
Shipping	Ice bag