

Anti-OLr1 Rabbit Polyclonal Antibody For Wb, Elisa - P78380

Item Number: CM0014440



Introduction

Detect oxidized LDL receptor 1 (OLR1/LOX-1) with this rabbit polyclonal antibody validated for WB and ELISA. Suitable for atherosclerosis and cardiovascular disease research. Cross-reacts with mouse OLR1.

[Learn More](#)

Parameter	Description
Product Name	OLR1 Rabbit pAb
Remarks/Alias	LOX1; LOXIN; SLOX1; CLEC8A; SCARE1; OLR1
Species	Human
GeneID (Human)	4973
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human OLR1 (NP_002534.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P78380
Protein Weight	31kDa
Shipping	Ice bag