

Anti-Oma1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96E52

Item Number: CM0019428

Introduction

Rabbit polyclonal antibody against human OMA1, a mitochondrial metalloprotease involved in quality control and apoptosis. Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for mitochondrial dynamics and stress response studies.

[Learn More](#)

Parameter	Value
Product Name	OMA1 Rabbit pAb
Remarks/Alias	DAB1; MPRP1; MPRP-1; YKR087C; ZMPOMA1; peptidase; 2010001O09Rik; OMA1
Species	Human
GeneID (Human)	115209
GeneID	115209
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 421-520 of human OMA1 (NP_660286.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96E52
Protein Weight	60kDa
Shipping	Ice bag