

Anti-Oma1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96E52

Item Number: CM0015988

Introduction

Rabbit polyclonal antibody targeting OMA1, a mitochondrial metalloprotease essential for quality control and stress response. Validated for WB and ELISA, cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	OMA1 Rabbit pAb
Remarks/Alias	DAB1; MPRP1; MPRP-1; YKR087C; ZMPOMA1; peptidase; 2010001009Rik; OMA1
Species	Human
GeneID (Human)	115209
GeneID	115209
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96E52
Protein Weight	60kDa
Shipping	Ice bag