

Anti-Oncomodulin-2 Polyclonal Antibody For Wb, Elisa - P0Ce71

Item Number: CM0019740

Introduction

Rabbit polyclonal antibody against human oncomodulin-2 (OM, OCM, OCM2). Validated for western blot and ELISA, cross-reacts with rat. Ideal for calcium-binding protein research.

[Learn More](#)

Parameter	Value
Product Name	OCM2 Rabbit pAb
Remarks/Alias	OM; OCM; OCM2
Species	Human
GeneID (Human)	4951
GeneID	4951
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-109 of human OCM2 (NP_006179.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P0CE71
Protein Weight	12kDa
Shipping	Ice bag