

Anti-Oncostatin M Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P13725

Item Number: CM0010628



Introduction

Recombinant rabbit monoclonal antibody against human Oncostatin M (OSM) validated for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studies in inflammation, cancer, and cytokine signaling.

[Learn More](#)

Parameter	Value
Product Name	Oncostatin M Rabbit mAb
Remarks/Alias	OSM; oncostatin-M; Oncostatin M
Species	Human
GeneID (Human)	5008
GeneID	5008
Immunogen	Synthetic peptide: a sequence within amino acids 1-100 of human Oncostatin M (P13725)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13725
Protein Weight	28kDa
Shipping	Ice bag