

Anti-Onecut1 Polyclonal Antibody For Wb - Q9Ubc0

Item Number: CM0015350

Introduction

Rabbit polyclonal antibody against human ONECUT1 (HNF6), suitable for WB and ELISA. Validated in human, mouse, and rat. Relevant for liver gene transcription studies.

[Learn More](#)

Parameter	Value
Product Name	ONECUT1 Rabbit pAb
Remarks/Alias	HNF6; HNF-6; HNF6A; ONECUT1
Species	Human
GeneID (Human)	3175
GeneID	3175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 167-292 of human ONECUT1 (NP_004489.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBC0
Protein Weight	51kDa
Shipping	Ice bag