


Anti-Opa1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, If-P, Elisa - O60313

Item Number: CM0009704



Introduction

OPA1 rabbit monoclonal antibody validated for WB, IP, IF/ICC, IF-P, and ELISA. Human OPA1 immunogen, cross-reacts with human, mouse, rat. For mitochondrial fusion, cristae remodeling, and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	OPA1 Rabbit mAb
Remarks/Alias	NPG; NTG; MGM1; BERHS; largeG; MTDPS14
Species	Human
GeneID (Human)	4976
GeneID	4976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-573 of human OPA1 (NP_056375.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60313
Protein Weight	108-118 kDa
Shipping	Ice bag