


Anti-Opa1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - O60313

Item Number: CM0025167



Introduction

Anti-OPA1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IHC-P, ELISA. Detects human, mouse, rat OPA1 (~112 kDa), a mitochondrial dynamin-like GTPase essential for inner membrane fusion and cristae remodeling.

[Learn More](#)

Parameter	Value
Product Name	OPA1 Rabbit pAb
Remarks/Alias	NPG; NTG; MGM1; BERHS; largeG; MTDPS14; OPA1
Species	Human
GeneID (human)	4976
GeneID	4976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-600 of human OPA1 (NP_056375.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60313
Protein Weight	112kDa
Shipping	Ice bag