

Anti-Opcml Polyclonal Antibody For Wb, Elisa - Q14982

Item Number: CM0030945



Introduction

High-quality rabbit polyclonal antibody targeting OPCML (opioid-binding protein/cell adhesion molecule, IgLON family member 1). Suitable for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat. Expected molecular weight: 38kDa. Ideal for neuroscience and cell adhesion research - Q14982.

[Learn More](#)

Parameter	Value
Product Name	OPCML Rabbit pAb
Remarks/Alias	OPCM; OBCAM; IGLON1; OPCML
Species	Human
GeneID	4978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 246-345 of human OPCML (NP_002536.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14982
Protein Weight	38kDa
Shipping	Ice bag