

Anti-Oplah Rabbit Polyclonal Antibody For Wb And Elisa - O14841

Item Number: CM0014790

Introduction

Rabbit polyclonal antibody against human OPLAH (5-oxoprolinase), validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for glutamate metabolism and glutathione degradation research.

[Learn More](#)

Parameter	Value
Product Name	OPLAH Rabbit pAb
Remarks/Alias	OPLA; OPLAHD; 5-Opase; OPLAH
Species	Human
GeneID (Human)	26873
GeneID	26873
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1119-1288 of human OPLAH (NP_060040.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14841
Protein Weight	137kDa
Shipping	Ice bag