

Anti-Opn5 Rabbit Polyclonal Antibody For Wb, Elisa - Q6U736

Item Number: CM0019964

Introduction

CamelBio's OPN5 rabbit polyclonal antibody, validated for WB and ELISA, cross-reacts with Human, Mouse, Rat. Ideal for studying opsin-5 photoreceptor function in circadian rhythms and retinal development.

[Learn More](#)

Parameter	Value
Product Name	OPN5 Rabbit pAb
Remarks/Alias	PGR12; GPR136; GRP136; TMEM13; OPN5
Species	Human
Gene ID (Human)	221391
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human OPN5 (NP_859528.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6U736
Protein Weight	40kDa
Shipping	Ice bag