

Anti-Oprd1 Rabbit Monoclonal Antibody For Wb, Elisa - P41143

Item Number: CM0004571



Introduction

High-quality rabbit monoclonal antibody specific to human OPRD1 (delta opioid receptor). Suitable for Western blot and ELISA applications, cross-reacts with mouse and rat. Essential for studying pain perception and opioid analgesia.

[Learn More](#)

Parameter	Value
Product Name	OPRD1 Rabbit mAb
Remarks/Alias	DOP; DOR; DOR1; OPRD
Species	Human
GeneID (Human)	4985
GeneID	4985
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human OPRD1 (P41143).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41143
Protein Weight	40kDa
Shipping	Ice bag