

Anti-Oprr1 Rabbit Polyclonal Antibody For Wb, Elisa - P41143

Item Number: CM0030049



Introduction

High-quality rabbit polyclonal antibody targeting human OPRD1 (delta-opioid receptor), validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for neuroscience and pain research applications.

[Learn More](#)

Parameter	Value
Product Name	OPRD1 Rabbit pAb
Remarks/Alias	DOP; DOR; DOR1; OPRD; OPRD1
Species	Human
GeneID (Human)	4985
GeneID	4985
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-372 of human OPRD1 (NP_000902.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41143
Protein Weight	40kDa
Shipping	Ice bag