

Anti-Oprk1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P41145

Item Number: CM0014842



Introduction

Rabbit polyclonal antibody targeting human OPRK1 (kappa opioid receptor), suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for pain and neuroscience studies.

[Learn More](#)

Parameter	Specification
Product Name	OPRK1 Rabbit pAb
Remarks/Alias	KOP; KOR; KOR1; OPRK; KOR-1; K-OR-1; OPRK1
Species	Human
GeneID	4986
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 280-380 of human OPRK1 (NP_000903.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41145
Protein Weight	43kDa
Shipping	Ice bag